

Using a Modified Signal Detection Framework in Obsessive-Compulsive Disorder (OCD): The OCD/UCM Bias Index

Richard Fiene PhD

Penn State Edna Bennett Pierce Prevention Research Center

August 2026

Abstract

This paper introduces a novel decision-making metric, the **OCD/UCM Bias Index**, designed to mathematically operationalize the "evidence gap" in Obsessive-Compulsive Disorder (OCD). While traditional models of OCD focus qualitatively on the intolerance of uncertainty, this framework utilizes a modified Signal Detection Theory (SDT) paradigm to isolate subjective response bias from baseline perceptual sensitivity. By evaluating an individual's decision-making thresholds under varying levels of stimulus ambiguity, the OCD/UCM Bias Index quantifies the exact volume of objective safety data required to counteract a subjective feeling of threat. By integrating both the Hit Rate (H) and the False Alarm Rate (F), this index mathematically captures the threat-avoidance "pathological lock" characteristic of OCD while eliminating diagnostic blind spots inherent to single-rate tracking.

Introduction

Obsessive-Compulsive Disorder (OCD) is characterized clinically by persistent, intrusive thoughts (obsessions) and repetitive physical or mental acts (compulsions) designed to mitigate distress. At the core of OCD's cognitive pathology lies a profound intolerance of uncertainty and an asymmetric "evidence gap"—a subjective state wherein an individual requires an extraordinary, often unattainable volume of objective data to neutralize an internal feeling of risk or error.^{1,2}

In classical cognitive-behavioral models, this phenomenon is described qualitatively. However, to advance the empirical evaluation of OCD treatments, we must mathematically operationalize how individuals navigate threat and safety under ambiguous conditions. This paper proposes the Uncertainty-Certainty Matrix (UCM)(See Table 1) as a conceptual framework³, bridged with the mathematical rigor of Signal Detection Theory (SDT), to cleanly quantify response biases in clinical populations.

The Threat-Detection Paradigm

To isolate and measure these cognitive biases, we propose a computerized laboratory task where participants are presented with a series of ambiguous visual stimuli (e.g., a visual gauge representing an environmental asset, such as an oven temperature dial). Across multiple trials, the signal-to-noise ratio of the gauge is manipulated to introduce varying levels of ambiguity.

To align this framework with standard psychological literature, the trial parameters are defined as follows:

- **Target Signal (Danger):** The gauge is objectively in an unsafe or faulty state (S_+).
- **Noise State (Safety):** The gauge is objectively in a safe or correct state (S_-).
- **Subjective Response (+):** The participant declares the system "Unsafe" (Threat detected).
- **Subjective Response (-):** The participant declares the system "Safe" (Threat dismissed).

Table 1: Standardized Uncertainty-Certainty Matrix (UCM)

	Subjective Response (-)"Safe" (Threat Dismissed)	Subjective Response (+)"Unsafe" (Threat Detected)
Noise State (S_-) Objective Safety	True Negative (TN) <i>Certainty Axis (Correct Rejection)</i>	False Positive (FP) <i>Uncertainty Axis (False Alarm)</i>
Signal State (S_+) Objective Danger	False Negative (FN) <i>Uncertainty Axis (Miss)</i>	True Positive (TP) <i>Certainty Axis (Hit)</i>

Primary Signal Detection Metrics & Mathematical Formalization

Standard SDT separates an individual's ability to discriminate between safety and danger from their overall willingness to report a threat. The behavioral responses are converted into standard rates:

- Hit Rate (H):** The probability of correctly identifying danger when danger is present.

$$\text{Hit Rate } (H) = \frac{TP}{TP + FN}$$

- False Alarm Rate (F):** The probability of misidentifying safety as danger.

$$\text{False Alarm Rate } (F) = \frac{FP}{TN + FP}$$

From these rates, we derive baseline perceptual sensitivity (d') and response criterion (c) using standard normal distribution transforms (Z-scores):

$$d' = Z(H) - Z(F)$$
$$c = -0.5 \times [Z(H) + Z(F)]$$

The OCD/UCM Bias Index

While traditional SDT metrics (d' and c) are mathematically robust, they require non-linear statistical conversions that complicate real-time clinical monitoring or straightforward digital dashboard interfaces. To deliver a linear, bounded behavioral coefficient that preserves the exact directional movement of psychological bias, we introduce the **OCD/UCM Bias Index**:

$$\text{OCD/UCM Bias Index} = 1 - (H + F)$$

By utilizing both foundational rates, this mathematical formulation scales behavior cleanly between -1.0 and $+1.0$. Incorporating the Hit Rate (H) alongside the False Alarm Rate (F) ensures that the metric accounts for the entire behavioral profile, resolving the mathematical limitations of earlier safe-trial models.

Clinical and Behavioral Interpretation

The mathematical behavior of the OCD/UCM Bias Index directly corresponds to distinct operational profiles across three core psychological states:

Optimal Calibration (Healthy Cohorts)

In non-clinical populations, individuals maintain a balanced response criterion ($c \approx 0$). Their visual detection systems register threats effectively while maintaining a minimal rate of false alarms. Under these conditions, the sum of the Hit Rate and False Alarm Rate hovers around parity:

$$H + F \approx 1.0 \Rightarrow \text{OCD/UCM Bias Index} \approx 0.0$$

A score of 0.0 mathematically represents ideal calibration, where internal perception functions in balance with objective reality.

The Pathological Lock (OCD Cohorts)

In threat-avoidance disorders like OCD, individuals exhibit a hyper-liberal response bias to prevent the catastrophic cost of a "Miss" (False Negative). As the threshold for reporting danger drops, the client treats mild ambiguity as certain danger, driving both the Hit Rate (H) and False Alarm Rate (F) toward maximum values ($H \rightarrow 1.0, F \rightarrow 1.0$). This causes the sum of the rates to exceed 1.5, forcing the index deep into negative values:

$$H + F > 1.5 \Rightarrow \text{OCD/UCM Bias Index} \ll 0.0$$

This negative range mathematically defines the **Pathological Lock**. The patient is trapped in the False Positive cell of the matrix, executing rituals to neutralize an environmental safety state that their cognitive threshold cannot register.

Conservative Bias

Conversely, if an individual or system demands absolute cognitive certainty before declaring a threat (e.g., a criminal jury threshold requiring proof "beyond a reasonable doubt"), both H and F contract toward zero. Under this hyper-conservative strategy:

$$H + F < 0.5 \Rightarrow \text{OCD/UCM Bias Index} \gg 0.0$$

This produces a positive index value approaching +1.0, capturing a strict refusal to trigger an alert without overwhelming evidentiary support.

Mathematical Soundness Assessment

This following section offers a practical and mathematical bridge between cognitive-behavioral theories of OCD and quantitative decision theory.

By introducing the **OCD/UCM Bias Index** ($\text{Bias Index} = 1 - [H + F]$), the paper presents a linear, non-parametric metric that operationalizes response bias without requiring normal inverse cumulative distribution transforms.

$$\text{OCD/UCM Bias Index} = 1 - (H + F)$$

Where:

- $H = \text{Hit Rate} = \frac{TP}{TP+FN}$
- $F = \text{False Alarm Rate} = \frac{FP}{TN+FP}$
- **Strict Boundedness:** Because $H, F \in [0,1]$, their sum $H + F \in [0,2]$. Consequently, $1 - (H + F)$ is strictly bounded in $[-1.0, +1.0]$.
- **Neutral Invariance** ($c = 0 \Leftrightarrow \text{Bias Index} = 0$): Under Gaussian Signal Detection Theory (SDT) assumptions:

$$H = \Phi\left(\frac{d'}{2} - c\right), \quad F = \Phi\left(-\frac{d'}{2} - c\right)$$

When an individual exhibits an unbiased, neutral response criterion ($c = 0$), $H = \Phi(d'/2)$ and $F = \Phi(-d'/2) = 1 - \Phi(d'/2)$. Thus:

$$H + F = \Phi\left(\frac{d'}{2}\right) + 1 - \Phi\left(\frac{d'}{2}\right) = 1.0 \Rightarrow 1 - (H + F) = 0.0$$

This property guarantees that **neutral decision-making always evaluates to 0.0, regardless of perceptual sensitivity (d')**.

- **Directional Symmetry:**
 - **Negative Values** (< 0.0): $H + F > 1.0$, signaling a hyper-liberal threat detection threshold (typical of the "Pathological Lock" in OCD).
 - **Positive Values** (> 0.0): $H + F < 1.0$, signaling a conservative/skeptical threshold.
- **Edge-Case Robustness:** Classical SDT response criterion ($c = -0.5[Z(H) + Z(F)]$) relies on probit transforms ($Z = \Phi^{-1}$), which become undefined ($\pm\infty$) when H or F equals 0.0 or 1.0 (requiring extreme-value correction log-edging like $\frac{1}{2N}$). The OCD/UCM Bias Index is continuously defined across the entire probability space without requiring artificial adjustments.

Mathematical Nuances & Limitations

Sensitivity Compression Effect: While Bias Index = 0.0 strictly aligns with $c = 0$, for non-zero criteria ($c \neq 0$), the scale magnitude of $1 - (H + F)$ is non-linearly attenuated by the individual's sensitivity (d'). Higher sensitivity shrinks $1 - (H + F)$ toward zero because both H and F approach 1.0 and 0.0 respectively.

While this does not impair diagnostic classification, researchers comparing raw Bias Index values across cohorts with vastly different perceptual discrimination capabilities should note that $1 - (H + F)$ acts as a probability-scale proxy for c rather than a latent normal Z-score.

Utility & Clinical Significance Table


Feature	Standard SDT (c)	Single-Rate Tracking (F or H alone)	OCD/UCM Bias Index
Mathematical Complexity	Non-linear / Inverse Normal (Z)	Linear	Linear / Basic Arithmetic
Edge-Case Stability	Fails at 0.0 and 1.0	Stable	Completely Stable
Sensitivity Independence	Fully Independent	Dependent	Exact at Neutral ($c = 0$)
Real-Time Dashboard Suitability	Low (Requires stats library)	High	High ($O(1)$ computation)
Interpretability for Clinicians	Abstract Z -units ($-\infty$ to $+\infty$)	Unidirectional Probability	Bounded Scale (-1.0 to $+1.0$)

Core Utility Advantages

- Digital Application Integration:** Its $O(1)$ arithmetic nature makes it computationally lightweight, making it suitable for real-time clinical dashboards, mobile assessments, and web applications.
- Operationalizing the "Evidence Gap":** Provides an objective, task-based behavioral metric that directly measures how much objective safety data (S_-) an OCD patient misclassifies as threat (S_+), complementing subjective self-report questionnaires.

Final Assessment

The OCD/UCM Bias Index represents a **sound, practical, and highly accessible improvement** over complex parametric modeling for clinical applications. It preserves the theoretical core of Signal Detection Theory while eliminating mathematical singularities and computational overhead.



References

1. Belloch, A., Fornés, G., Carrasco, A., López-Solá, C., Alonso, P., & Menchón, J. M. (2016). Incompleteness and not just right experiences in the explanation of Obsessive-Compulsive Disorder. *Psychiatry Research*, 236, 1-8. <https://doi.org/10.1016/j.psychres.2016.01.012>
2. Carroll, I. C., Planalp, E. M., Van Hulle, C. A., & Goldsmith, H. H. (2019). Peer victimization and selective attention in adolescence: Evidence from a monozygotic twin difference design. *Journal of Abnormal Child Psychology*, 47(8), 1303-1313. <https://doi.org/10.1007/s10802-019-00516-7>
3. Fiene, R. (2025). Uncertainty-Certainty Matrix. *Knowledge Journal*, 5(2), 8. <https://doi.org/10.3390/knowledge5020008>

Updated OCD/UCM Bias Index Calculator v4

Step 1: Trial Data Entry

Metric	Count
True Negative (TN)	
False Positive (FP)	
False Negative (FN)	
True Positive (TP)	

Step 2: SDT Metrics & Bias Index

Metric	Formula / Rate
Hit Rate (Raw, H)	
False Alarm Rate (Raw, F)	
Hit Rate (Corrected)	
FA Rate (Corrected)	
Sensitivity (d')	
Criterion (c)	
OCD/UCM Bias Index	

Step 3: Clinical Diagnostic Profile

Diagnostic Category	Profile Readout
Clinical State	

OCD/UCM Bias Index Scale

Pathological Lock = -.4+

Conservative Bias = +.4+


Optimal Calibration = -.39 -- +.39

See 15 simulated cases generated below from the OCD/UCM Bias Index Calculator.

Fifteen (15) simulated examples calculating the OCD/UCM Bias Index across the three diagnostic profiles using the OCD/UCM Bias Index: $1 - (H + F)$.


OCD/UCM Bias Index: Trial Data & Calculations



ID	(TN)	(FP)	(FN)	(TP)	Hit Rate (H)	FA Rate (F)	Bias Index	Diagnostic Profile
#1	5	45	2	48	0.960	0.900	-0.860	Pathological Lock
#2	10	40	5	45	0.900	0.800	-0.700	Pathological Lock
#3	15	35	1	49	0.980	0.700	-0.680	Pathological Lock
#4	8	42	10	40	0.800	0.840	-0.640	Pathological Lock
#5	2	48	0	50	0.990	0.960	-0.950	Pathological Lock
#6	48	2	40	10	0.200	0.040	+0.760	Conservative Bias
#7	45	5	45	5	0.100	0.100	+0.800	Conservative Bias
#8	50	0	38	12	0.240	0.010	+0.750	Conservative Bias
#9	42	8	35	15	0.300	0.160	+0.540	Conservative Bias
#10	49	1	48	2	0.040	0.020	+0.940	Conservative Bias
#11	45	5	5	45	0.900	0.100	0.000	Optimal Calibration
#12	40	10	10	40	0.800	0.200	0.000	Optimal Calibration
#13	48	2	8	42	0.840	0.040	+0.120	Optimal Calibration
#14	45	5	15	35	0.700	0.100	+0.200	Optimal Calibration
#15	45	5	20	30	0.600	0.100	+0.300	Optimal Calibration




The above graph provides the results from 500 Bias Index cases in just providing the OCD/UCM Bias Index Coefficients

Three Examples of OCD/UCM Bias Indices: Positive, Negative, Neutral





Clinical ERP (Exposure & Response Prevention) Daily Progress Dashboard



The above graph presents how all the data from the OCD/UCM Bias Index and Application could potentially be displayed for clinical use in an ERP Treatment Plan.